

PCR Primers

Gene	Primer sequences	Species *	#cycles	Product size
Insulin	TCCTACCCCTGCTGGCCCTGC AGTTGCAGTAGTTCTCCAG	m	18	305
nkx6.1	TCTTCTGGCCTGGGGTGATG GTGCTTCTTTCTCCACTTGGTCC	r	27	277
nkx2.2	CATCTTGGACCTTCCGGACAC GGACTTGGAGCTCGAGTCTTG	m	27	267
pdx-1	ACCACCTTCCAGCTCAGCTCC CGGGTGTAGGCAGTACGGGTC	m	30	213
Pax4	ACCCTGTGACATTTACGGAG CGATTGATAGAGGACACACTGG	m	30	260
Pax6	TCACAGCGGAGTGAATCAGCT TACTCACACAACCGTTGGATACT	m	27	204
BETA2/NeuroD	CTTGGCCAAGAATACTATCTGG GGAGTAGGGATGCACCGGAA	m	25	228
glucokinase	TGAGCCCAGAGAAGAAAGCT CTGGAACCTGCCAGGAT	m	30	479
GLUT2	CAAGATCACCGAACCTTGG ATCCGCCCACTGCAAAGCT	m	30	311
SUR1	CAGTCATGGAGAGGAAAGCCCC GGATGATGCGATTGAGGAGGCT	m	27	493
α -tubulin	GCGTGAGTGTATCTCTATCCAC GGTAGGTTCCGGTGCGTACTT	m	18	247
G6PD	ACCTGCAGAGCTCCAATCAAC CACACAAGCAATGTTGTCTCG	m	30	172
SV40 T antigen	ATCCTGGTGTGATGCAATG CAGAGAGGAATCTTTCAGC		30	272
Isl-1	AGCAAGAACGACTTCGTGATG GACTGAGAGGGTCTCCAGCTC	r/h	27	187
HNF3 β	CGCTCGGGACCCCAAGACATACC GGCGAGCGGGGCACCTTGAGAAA	m	30	240
HNF6	TCCCGCGCACCCGCGTTCCATG CGCATAAGTGTGAAACTACCGC	m	30	193

* m=mouse; r=rat; h=hamster

TOUCHDOWN RT-PCR PROTOCOL

This protocol uses the **Titan™ One Tube RT-PCR System** from Boehringer Mannheim & a Perkin Elmer 9700 PCR System. Primers are at 10 μ M.

In two separate tubes, prepare the following:

1 μ l top primer
1 μ l bottom primer
1 μ l α -tubulin top (internal control)
1 μ l α -tubulin bottom (internal control)
1 μ l 10 mM dATP (USB)
0.5 μ l 10 mM dCTP (USB)
1 μ l 10 mM dGTP (USB)
1 μ l 10 mM dTTP (USB)
0.25 μ l [α -³²P] dCTP
2.5 μ l DTT solution (kit)
0.25 μ l Rnase Inhibitors (Promega)
50 ng RNA template
water to 25 μ l V₁

14 μ l water
10 μ l 5X RT-PCR buffer (kit)
1 μ l enzyme mix (kit)
water to 25 μ l V₂

Pre-warm PCR machine to 50°C. Combine tubes 1 & 2 in a 0.2 μ l PCR reaction tube and cycle as follows:

Step 1	50° x 30 min; 94° x 2 min	(1 cycle)
Step 2	94° x 1 min; 63° x 1 min; 72° x 1 min	(2 cycles)
Step 3	94° x 1 min; 62° x 1 min; 72° x 1 min	(2 cycles)
Step 4	94° x 1 min; 61° x 1 min; 72° x 1 min	(2 cycles)
Step 5	94° x 1 min; 60° x 1 min; 72° x 1 min	(2 cycles)
Step 6	94° x 1 min; 59° x 1 min; 72° x 1 min	(34 cycles)
Step 7	72° x 7 min; 4° hold	

Remove samples at appropriate time points (i.e., cycles) and add to an equal volume of formamide loading dye.

Load samples in a 6% sequencing gel and run until first dye band runs off bottom of gel. Dry gel and expose to film overnight.